

MIN Faculty Department of Informatics



Grasping with the Shadow Dexterous Hand integrated on the mobile robot PR2

Yannick Jonetzko



University of Hamburg Faculty of Mathematics, Informatics and Natural Sciences Department of Informatics

Technical Aspects of Multimodal Systems

11. July 2017



Outline

2	outime				
	Introduction	Shadow Hand	BioTac	Stable grasp	Integration with the PR2
	1. Introducti	on			
	2. Shadow H	land			
	3. BioTac				
	Raw D	ata			
	Filter				
	Norma	lization			
	Point o	of Contact			
	Calibra	ation			
	4. Stable gra	asp			
	Contro	oller			

5. Integration with the PR2



BioTac

Stable g

Integration with the PR2

- Service robots for daily tasks
 - \rightarrow Picking up objects
- Anthropomorphic gripper in the human environment
- One gripper for "all" kinds of objects





BioTac

Stable grasp

Integration with the PR2

- Project 2015: Delivery service
- Student Assistant: Different tasks of grasping
 - Grasping action server for Robotiq gripper
 - More than just a simple gripper



robotiq s model gripper



Designed to imitate and reproduce dexterity of a human hand

- Size and shape
- Kinematics

Different models:

- Left, right hand
- Motor, air muscle controlled
- Amount of fingers
- BioTac sensors



Shadow Dexterous Hand

Shadow Hand



Shadow Dexterous Hand kinematics[src13]

Shadow Hand

Introduction	Shadow Hand	BioTac	Stable grasp	Integration with the PR2
M.L.				
video				I

https://www.youtube.com/watch?v=9ubXFMklEe8

Human sense:

- Cutaneous sensing: Pressure, temperature, pain
- ► Kinesthetic sensing: Position, movement, equilibrium

Measuring cutaneous sensing data with tactile sensors

Different types of tactile sensors:

- Force/Torque: Pressure and contact point
- Dynamic: Sliding, making and breaking contact
- Thermal: Temperature and thermal conductivity

BioTac cross section[st:15]

 Material identification:
 A. Gómez Eguíluz, I. Rañó, S.A. Coleman, T.M. McGinnity [ERCM16]

BioTac

- Radius of curvature: Nicholas Wettels, Gerald E. Loeb [WL11]
- Handover

A. Gómez Eguíluz, I. Rañó, S.A. Coleman, T.M. McGinnity [ERCM16]

Point of contact
 V. Ciovanu, D. Popescu, A. Petrescu [CPP14]

- ► Electrode voltages (E1 E19) Measure voltage
- ► DC Pressure (Pdc) Absolute fluid pressure
- ► AC Pressure (Pac) Dynamic fluid pressure (vibrations)
- DC Temperature (Tdc) Device temperature
- ► AC Temperature (Tac) Thermal flux

All incoming values from 0 to 4095

Introduction

Goal

Get information which are useful for grasping.

E.g.: Pressure (force), the point of contact

To reach this goal the data has to be prepared:

- filtered
- normalized
- calibrated (comparable, units)

Electrode Voltages (E1 - E19) - Impedance

DC Pressure (Pdc)

Video

Temperature (Tdc - Tac)

Introduction

Shadow Hand

BioTac

Stable gr

Video

Problems

Introduction

adow Hand

BioTac

Stable gras

- Different amount of liquid in every fingertip
- Temperature inflates volume
- Skin wear inflates volume
- $\rightarrow\,$ frequent calibration is necessary

Introduction Shadow Hand BioTac Stable grasp Integration with the PR2

Noisefiltering with the Exponentially Weighted Moving Average (EMA)

$$EMA_n = \alpha \cdot x_n + (1 - \alpha) \cdot EMA_{n-1}$$

With:

- \blacktriangleright EMA₀ = 0
- ► x_n: New value
- Advantage: Reduce noise, smoothing the curves
- Disadvantage: Delay of values with low alpha

All values should be comparable to each other.

Naive approach: Set all values to 0 at the beginning.

 \rightarrow Problem: The values are drifting.

Better approach: Detect contacts and zeroing values when there is no contact.

BioTac

Stable

Integration with the PR2

- Assumption: Starting the normalization in a rest state
- Conditions for no contact:
 - Mean of Pdc values has to be under threshold (5).
 - Difference between min and max value of electrodes has to be under threshold (10).
- Prevent zeroing with very small contacts:
 - $\rightarrow\,$ Safe values every 0.5 seconds and mean over 5 values.

	BioTac	

Possible solutions:

- Improve thresholds
- Consider finger efforts
- External service when there is no contact (e.g. after grasp)

Point of Contact

		BioTac	
Video			
	bytemat the second se	Apply cam/intege_rev 2	
	Ad Declars in		

Point of Contact

Introduction

hadow Hand

BioTac

Stable gras

Idea:

- 1. Get the average of the electrode position, weighted with the impedance value
- $\rightarrow\,$ position on the core of the BioTac
- 2. Map position to surface and distinguish between:
 - Cylindrical part
 - Spherical part

[CHLL13]

Compute the PoC by a weighted average of the position of the electrodes:

BioTac

$$\langle x_c, y_c, z_c \rangle = \frac{\sum_{i=1}^{19} (|e_{i^*}|^n \langle x_i, y_i, z_i \rangle)}{\sum_{i=1}^{19} (|e_{i^*}|^n)}$$

*e_i**: Normalized electrode values *n*: Impect of the single electrode

 $\rightarrow\,$ Turns out, considering just the negative values produces a way better result.

Point of Contact Map to surface

		BioTac	
Cylindrica	l part		
$\langle x_{c'}, y_{c'}, z_{c'} \rangle$	$\langle x_c, \frac{r*y_c}{\sqrt{y_c^2+z}} \rangle = \langle x_c, \frac{r*y_c}{\sqrt{y_c^2+z}} \rangle$	$\frac{1}{\sqrt{2}}, \frac{r*z_c}{\sqrt{y_c^2+z_c^2}}$	

Spherical part

$$\langle x_{c'}, y_{c'}, z_{c'} \rangle = \langle x_c, y_c, z_c \rangle \frac{r}{\sqrt{x_c^2 + y_c^2 + z_c^2}}$$

[CHLL13]

 Introduction
 Shadow Hand
 BioTac
 Stable grasp
 Integration with the I

 Positioning similar to point of contact, without weighted average

Color:

- Red: Small impedance
- Blue: High impedance

Current state

For a stable grasp the calibrated pressure is the most interesting part:

- 1. Calibrate fingers among each other
 - \rightarrow Press finger at the same contact point against each other.
- 2. Calibrate with real force units
 - → External way to measure force is necessary. It might work with effort values from Shadow Hand fingers.

http://wiki.ros.org/ros_control

Stable grasp

Integration with the PR2

- Effort controller
- Position controller
- Velocity controller
- Trajectory controller
- Mixed controller

Goal

Integrate tactile feedback to hold an effort.

Possible solution:

- 1. Close fingers with effort and velocity controlled motion.
- 2. Stop motion when contact is detected.
- **3**. Hold effort within a threshold when two contact points with opposite force directions are detected.

Problems: Moving hand around affects pressure feedback.

Integration with the PR2

Shadow Hand

BioTac

Stable gras

Integration with the PR2

Current state:

- Simulation works
- Shadow Hand is controlled with base station
- Shadow Hand and PR2 are not tested together

References

Integration with the PR2 [CHLL13] Jeremy A. Fishel Chia-Hsien Lin and Gerald E. Loeb. Estimating point of contact, force and torque in a biomimetic tactile sensor with deformable skin SynTouch LLC, 2013. [CPP14] Vlad Ciobanu, Decebal Popescu, and Adrian Petrescu. Point of contact location and normal force estimation using biomimetical tactile sensors In 2014 Eighth International Conference on Complex, Intelligent and Software Intensive Systems. IEEE, jul 2014. [Dev04] Developments in dextrous hands for advanced robotic applications. Automation Congress, Jun 2004. [ERCM16] A. Gomez Eguiluz, I. Rano, S.A. Coleman, and T.M. McGinnity. A multi-modal approach to continuous material identification through

tactile sensing.

In 2016 IEEE/RSJ International Conference on Intelligent Robots and Systems (IROS). IEEE, oct 2016.

References (cont.)

Introduction

- [NP13] Kien-Cuong Nguyen and Veronique Perdereau. Fingertip force control based on max torque adjustment for dexterous manipulation of an anthropomorphic hand. In 2013 IEEE/RSJ International Conference on Intelligent Robots and Systems. IEEE, nov 2013.
 [src13] Shadow Dexterous Hand Technical Specification, 1 2013.
- [st:15] BioTac Product Manual, 3 2015.
- [WL11] Nicholas Wettels and Gerald E. Loeb.
 Haptic feature extraction from a biomimetic tactile sensor: Force, contact location and curvature.
 In 2011 IEEE International Conference on Robotics and Biomimetics.
 IEEE, dec 2011.